




Tel: 023 95 49 86

Mail: contact@ghl-medical.com / comercial@ghl-medical.com

Adresse: SARL GHL MEDICAL N°53 SI EL HAOUAS AIN TAYA ALGER



- **BIPOLAR PROSTHESIS
(FENESTRATED)**

[PRÓTESIS BIPOLAR (FENESTRADO)]

SIZE	37 MM TO 55 MM (DIFF. 1 MM)
MATERIAL	SS



- **BIPOLAR PROSTHESIS
(NON-FENESTRATED)**

[PRÓTESIS BIPOLAR (NO FENESTRADO)]

SIZE	37 MM TO 55 MM (DIFF. 1 MM)
MATERIAL	SS



● AUSTIN MOORE HIP PROSTHESIS

[PRÓTESIS DE CADERA AUSTIN MOORE]

SIZE	37 MM TO 55 MM (DIFF. 1 MM)
MATERIAL	SS



● AUSTIN MOORE HIP PROSTHESIS (NARROW STEM)

[PRÓTESIS DE CADERA AUSTIN MOORE (VAPOR ESTRECHO)]

SIZE	37 MM TO 55 MM (DIFF. 1 MM)
MATERIAL	SS



- **AUSTIN MOORE HIP PROSTHESIS (LONG STEM)**

[PRÓTESIS DE CADERA AUSTIN MOORE (TALLO LARGO)]

SIZE	37 MM TO 55 MM (DIFF. 1 MM)
MATERIAL	SS



- **THOMPSON HIP PROSTHESIS**

[PRÓTESIS DE CADERA THOMPSON]

SIZE	37 MM TO 55 MM (DIFF. 1 MM)
MATERIAL	SS



● THOMPSON HIP PROSTHESIS STANDARD

(LONG STEM)

[PRÓTESIS DE CADERA THOMPSON ESTÁNDAR
(TALLO LARGO)]

SIZE	37 MM TO 55 MM (DIFF. 1 MM)
MATERIAL	SS



● SHOULDER PROSTHESIS

[PRÓTESIS DE HOMBRO]

SIZE	31 MM TO 41 MM (DIFF. 2 MM)
MATERIAL	SS
DIA Ø	6.0, 0.7, 8.0 MM



FIX BIPOLAR

● CEMENTED STEM (POLISHED)

[VÁSTAGO CEMENTADO (PULIDO)]

MATERIAL	SIZE	LENGTH
SS 316L	NO - 8	115 MM
	NO - 9	130 MM
	NO - 10	140 MM
	NO - 11	145 MM
	NO - 12	150 MM
	NO - 13	155 MM
	NO - 14	160 MM
	NO - 15	165 MM

(12 / 14 TAPPER & NECK ANGLE 135°)
(SELF CENTERING)



● UNCEMENTED STEM (H. A. COATED)

[VÁSTAGO NO CEMENTADO (REVESTIDO H.A)]



MATERIAL	SIZE	LENGTH
TATINUM	NO - 8	115 MM
	NO - 9	130 MM
	NO - 10	140 MM
	NO - 11	145 MM
	NO - 12	150 MM
	NO - 13	155 MM
	NO - 14	160 MM
	NO - 15	165 MM

(12 / 14 TAPPER & NECK ANGLE 135°)
(SELF CENTERING)

FIX BIPOLAR

● UNCEMENTED STEM (POROUS COATED)

[VÁSTAGO NO CEMENTADO (REVESTIMIENTO POROSO)]

MATERIAL	SIZE	LENGTH
SS 316L	NO - 8	115 MM
	NO - 9	130 MM
	NO - 10	140 MM
	NO - 11	145 MM
	NO - 12	150 MM
	NO - 13	155 MM
	NO - 14	160 MM
	NO - 15	165 MM

(12 / 14 TAPPER & NECK ANGLE 135°)
(SELF CENTERING)




● MODULAR BIPOLAR CUP

[COPA BIPOLAR MODULAR]

SIZE	38 MM TO 56 MM
MATERIAL	SS 316L

(SELF CENTERING)



TOTAL HIP REPLACEMENT

● INNER HEAD (POLISHED)

[CABEZA INTERIOR (PULIDA)]

MATERIAL	DIA - Ø 28 MM (12 / 14 TAPPER)			
SS 316L	-3.5	0.0	+3.5	+7.0



● CENTRELIZER

[CENTRIZADOR]

● CEMENT RETRACTOR

[RETRACTOR DE CEMENTO]



TOTAL HIP REPLACEMENT

● UNCEMENTED ACETABULAR CUP WITH LINER UHMWPW 10° & 180° (H. A. COATED)

[COPA ACETABULAR NO CEMADA
CON LINER UHMWPW 10°Y 180°(H. A. RECUBIERTO)]

SIZE	44 MM, 46 MM, 48 MM, 50 MM, 52 MM, 54 MM, 56 MM, 58 MM, 60 MM, 62 MM, 64 MM
MATERIAL	TITANIUM (H.A. COATED)

(SELF CENTERING)



● UNCEMENTED ACETABULAR CUP WITH LINER UHMWPW 10° & 180° (POROUS COATED)

[COPA ACETABULAR NO CEMADA
CON LINER UHMWPW 10°Y 180°(REVESTIMIENTO POROSO)]



SIZE	44 MM, 46 MM, 48 MM, 50 MM, 52 MM, 54 MM, 56 MM, 58 MM, 60 MM, 62 MM, 64 MM
MATERIAL	TITANIUM (POROUS COATED)

(SELF CENTERING)

TOTAL HIP REPLACEMENT

● CEMENTED ACETABULAR CUP 10° & 180°

[COPA ACETABULAR CEMENTADA 10°Y 180°]

SIZE	44 MM, 46 MM, 48 MM, 50 MM, 52 MM, 54 MM, 56 MM, 58 MM, 60 MM, 62 MM, 64 MM
MATERIAL	UHMWPW

(SELF CENTERING)





● ACETABULAR CUP SCREW

[TORNILLO COPA ACETABULAR]

SIZE Ø6.5 MM	20 MM TO 50 MM (DIFF. 5 MM)
MATERIAL	TITANIUM



INSTRUMENT SET



TOTAL HIP REPLACEMENT





Contactez Nous

