



Tel: 023 95 49 86

Mail: contact@ghl-medical.com / comercial@ghl-medical.com

Adresse: SARL GHL MEDICAL N°53 SI EL HAOUAS AIN TAYA ALGER



Intramedullary Nails

● HUMERUS NAIL CANNULATED

[CLAVO DE HÚMERO CANULADO]

Ø6.5 MM - 3.4 MM Interlocking Bolt

Ø7.5 MM - 3.4 MM Interlocking Bolt

Distal Holes

3 Distal Holes

3 Proximal Holes

Proximal Holes

Ø6.5 MM - 3.4 MM Interlocking Bolt

Ø7.5 MM - 3.4 MM Interlocking Bolt

LENGTH	200, 220, 240, 260, 280, 300 MM
TYPE	SS / TT
DIA Ø	6.5, 7.3 MM

● TIBIA NAIL CANNULATED

[CLAVO TIBIA CANULADO]

Ø8.0 MM - 3.9 MM Interlocking Bolt

Ø9.0 MM to 11.0 MM - 4.9 MM Interlocking Bolt

Distal Holes

3 Distal Holes

3 Proximal Holes

Proximal Holes

Ø8.0 MM - 3.9 MM Interlocking Bolt

Ø9.0 MM to 11.0 MM - 4.9 MM Interlocking Bolt

LENGTH	280, 300, 320, 340, 360, 380 MM
TYPE	SS / TT
DIA Ø	8.0, 9.0, 10.0, 11.0 MM

● EXPERT TIBIA NAIL CANNULATED

[EXPERT TIBIA UÑA CANULADA]

Ø8.0 MM - 3.9 MM Interlocking Bolt

Ø9.0 MM to 11.0 MM - 3.9 & 4.9 MM Interlocking Bolt

Distal Holes

4 Distal Holes

5 Proximal Holes

Proximal Holes

Ø8.0 MM - 3.9 MM Interlocking Bolt

Ø9.0 to 11.0 MM - 4.9 MM Interlocking Bolt

LENGTH	280, 300, 320, 340, 360, 380 MM
TYPE	SS / TT
DIA Ø	8.0, 9.0, 10.0, 11.0 MM



Intramedullary Nails

● FEMUR NAIL CANNULATED

[CLAVO FÉMUR CANULADO]

Ø9.0 MM to 12.0 MM - 4.9 MM Interlocking Bolt

Distal Holes

2 Distal Holes

2 Proximal Holes

Proximal Holes

Ø9.0 MM to 12.0 MM - 4.9 MM Interlocking Bolt

LENGTH	320, 340, 360, 380, 400, 420, 440 MM
TYPE	SS / TT
DIA Ø	9.0, 10.0, 11.0, 12.0 MM

● 130° & 135° SHORT PROXIMAL FEMUR NAIL CANNULATED

[PROXIMAL CORTO 130° Y 135°
CLAVO FÉMUR CANULADO]

Ø8.0 MM - 3.9 MM Interlocking Bolt

Ø9.0 MM to 11.0 MM - 3.9 & 4.9 MM Interlocking Bolt

Distal Holes

4 Distal Holes

2 Proximal Holes

Proximal Holes

Ø9.0 MM TO 12.0 MM - 6.4 MM CANNULATED
PFN HIP SELF TAPPING

Ø9.0 MM TO 12.0 MM - 8.0 MM CANNULATED PFN
FEMORAL NECK SCREW
SELF TAPPING

LENGTH	280, 300, 320, 340, 360, 380 MM
TYPE	SS / TT
DIA Ø	8.0, 9.0, 10.0, 11.0 MM

● 130° & 135° TROCHANTERIC PROXIMAL FEMUR NAIL CANNULATED

[PROXIMAL TROCANTERICO DE 130° Y 135°
CLAVO FÉMUR CANULADO]

Ø9.0 MM to 12.0 MM - 4.9 MM Interlocking Bolt

Distal Holes

2 Distal Holes

2 Proximal Holes

Proximal Holes

Ø9.0 MM TO 12.0 MM - 6.4 MM CANNULATED PFN HIP SELF TAPPING

Ø9.0 MM TO 12.0 MM - 8.0 MM CANNULATED PFN FEMORAL NECK SCREW SELF TAPPING

LENGTH	180 MM
TYPE	SS / TT
DIA Ø	9.0, 10.0, 11.0, 12.0 MM




Intramedullary Nails

● 130° SHORT PROXIMAL FEMUR NAIL ANTIROTATION - II CANNULATED

[CLAVO FÉMUR PROXIMAL CORTO 130°
ANTIROTACIÓN - II CANULADO]


Ø9.0 MM to 12.0 MM - 4.9 MM Interlocking Bolt



● 130° TROCHANTERIC PROXIMAL FEMUR NAIL ANTIROTATION - II CANNULATED

[FÉMUR PROXIMAL TROCANTERICO 130°
UÑA ANTIROTACIÓN - II CANULADA]

Ø9.0 MM to 12.0 MM - 4.9 MM Interlocking Bolt
Distal Holes



● 130° & 135° LONG PROXIMAL FEMUR NAIL CANNULATED

[FÉMUR PROXIMAL LARGO DE 130° Y 135°
UÑA CANULADA]

Ø9.0 MM to 12.0 MM - 4.9 MM Interlocking Bolt



Proximal Holes

Ø9.0 MM TO 12.0 MM - 6.4 MM CANNULATED PFN HIP SELF TAPPING
Ø9.0 MM TO 12.0 MM - 8.0 MM CANNULATED PFN FEMORAL NECK
SCREW SELF TAPPING

LENGTH	340, 360, 380, 400, 420, 440 MM
TYPE	SS / TT
DIA Ø	9.0, 10.0, 11.0, 12.0 MM
ANGLED	LEFT / RIGHT

● 130° LONG PROXIMAL FEMUR NAIL ANTIROTATION - II CANNULATED

[CLAVO FEMUR PROXIMAL LARGO 130°
ANTIROTACIÓN - II CANULADO]

Ø9.0 MM to 12.0 MM - 4.9 MM Interlocking Bolt



Proximal Holes

Ø9.0 MM to 12.0 MM - 12.0 MM CANNULATED PFNA
- II BLADE

LENGTH	340, 360, 380, 400, 420, 440 MM
TYPE	SS / TT
DIA Ø	9.0, 10.0, 11.0, 12.0 MM
ANGLED	LEFT / RIGHT

Intramedullary Nails

● SUPRA CONDYLAR DISTAL FEMUR NAIL CANNULATED

[FÉMUR DISTAL SUPRACONDILAR
UÑA CANULADA]



LENGTH	340, 360, 380, 400, 420, 440 MM
TYPE	SS / TT
DIA Ø	9.0, 10.0, 11.0, 12.0 MM
ANGLED	LEFT / RIGHT

● ELASTIC NAIL (TENS)

[UÑA ELÁSTICA (TENS)]



1.5 MM
2.0 MM
2.5 MM
3.0 MM
3.5 MM
4.0 MM
4.5 MM

LENGTH	340, 360, 380, 400, 420, 440 MM
TYPE	SS / TT
DIA Ø	9.0, 10.0, 11.0, 12.0 MM
ANGLED	LEFT / RIGHT



Intramedullary Nails

- **ENDER'S NAIL**

[UÑA DE ENDER]



LENGTH	200, 220, 240, 260, 280, 300, 320, 340, 360, 380, 400, 420, 440 MM
TYPE	SS
DIA Ø	2.0, 2.5, 3.0, 3.5, 4.0, 4.5 MM

- **12.0 MM CANNULATED PFNA - II BLADE**


[PFNA CANULADO DE 12,0 MM - HOJA II]



LENGTH	200, 220, 240, 260, 280, 300, 320, 340, 360, 380, 400, 420, 440 MM
TYPE	SS
DIA Ø	2.0, 2.5, 3.0, 3.5, 4.0, 4.5 MM

- **SQUARE NAIL FOR ULNA**

[CLAVO CUADRADO PARA CÚBITO]



LENGTH	200, 220, 240, 260, 280, 300 MM
TYPE	SS
DIA Ø	2.0, 2.5, 3.0, 3.5, 4.0 MM



Intramedullary Nails

● SQUARE NAIL FOR RADIUS

[CLAVO CUADRADO PARA RADIO]



LENGTH	200, 220, 240, 260, 280, 300 MM
TYPE	SS
DIA Ø	2.0, 2.5, 3.0, 3.5, 4.0 MM

● RUSH NAIL


[CLAVO RÁPIDO]



LENGTH	200, 220, 240, 260, 280, 300 MM
TYPE	SS
DIA Ø	2.0, 2.5, 3.0, 3.5, 4.0 MM

● K - NAIL

[K-CLAVO]



LENGTH	280, 300, 320, 340, 360, 380 MM
TYPE	SS / TT
DIA Ø	8.0, 9.0, 10.0, 11.0, 12.0, 13.0 MM



Intramedullary Nails

- **END CAP - TIBIA / FEMUR / PEN / SUPRA**

[TAPA DE EXTREMO - TIBIA / FEMUR / PEN / SUPRA]



LENGTH	0, 5, 10, 15 MM
TYPE	SS / TT

- **3.4 MM INTERLOCKING BOLT SELF TAPPING**

[TORNILLO DE ENCLAVAMIENTO AUTORROSCANTE DE 3,4 MM]



LENGTH	14 TO 50 MM (DIFF. 2 MM)
TYPE	SS / TT
DRILL BIT	2.7 MM

- **3.9 MM INTERLOCKING BOLT SELF TAPPING**

[TORNILLO DE ENCLAVAMIENTO AUTORROSCANTE DE 3,9 MM]



LENGTH	14 TO 80 MM (DIFF. 2 MM)
TYPE	SS / TT
DRILL BIT	3.2 MM



Intramedullary Nails

● 4.9 MM INTERLOCKING BOLT SELF TAPPING

[TORNILLO DE ENCLAVAMIENTO AUTORROSCANTE DE 4,9 MM]



LENGTH	18 TO 80 MM (DIFF. 5 MM)
TYPE	SS / TT
DRILL BIT	4.0 MM

● 6.4 MM CANNULATED PFN HIP PIN SELF TAPPING

[PERNO DE CADERA PFN CANULADO DE 6,4 MM AUTORROSCANTE]



LENGTH	40 TO 110 MM (DIFF. 5 MM)
TYPE	SS / TT
REAMER	6.4 MM
GUIDE WIRE	2.5 MM

● 8.0 MM CANNULATED PFN FEMORAL NECK SCREW SELF TAPPING

[TORNILLO DE CUELLO FEMORAL PFN CANULADO DE 8,0 MM AUTOPERFORANTES]




LENGTH	40 TO 110 MM (DIFF. 5 MM)
TYPE	SS / TT
REAMER	8.0 MM
GUIDE WIRE	2.5 MM



PFN II INSTRUMENT SET


   FOLLOW US



**PFN II
INSTRUMENT SET**

TENS NAIL INSTRUMENT

TENS NAIL **INSTRUMENT SYSTEM**




HUMERUS NAIL INSTRUMENT



   FOLLOW US

HUMERUS NAIL INSTRUMENT SYSTEM



CONTACT US



Contactez Nous

